

# How To Order the Infinite Socket™ System

## Step 1 - Register

- Register for an account or login at "[liminnovations.com/login](https://liminnovations.com/login)"

## Step 2 - Order Form

- After you login, click "Order" to pull up the online order form.
- Fill out the online Order Form in its entirety. Incomplete forms and offline order submissions will delay processing.
- Use the "measurement form," available for download under "Your Account" as a guide. *\*All measurements are required for order processing\**

## Step 3 - Order Method (Choose one option to accompany your online order)

### A. LIM Capture - Digital

Compatible Sockets:  
TF Pin-lock/lanyard.

1. Take Frontal and Sagittal view photos,
2. Photos are taken 8 feet away at level of umbrella.




3. Upload the images to the online order form when prompted
4. Submit online order.

### B. Scans

Compatible Sockets:  
TF/TT all suspensions.

1. Scan the residual limb.



2. Save the file in STL format.
3. Email scan to [scans@liminnovations.com](mailto:scans@liminnovations.com) along with patient name.
4. Submit online order.

### C. Limb Impression

Compatible Sockets:  
TF/TT all suspensions.

1. Take a fiberglass or plaster mould of the patients residual limb.
2. Mark bony landmarks and alignment lines



3. Mail mold to LIM Innovations®.  
*\*See address below\**
4. Submit online order.

For any questions, please contact customer service: (844) 888-8546

Follow LIM Innovations®



LIMinnovations.com  
[info@liminnovations.com](mailto:info@liminnovations.com)  
844.888.8546

424 9th Street  
San Francisco  
CA 94103

