

THE INFINITE SOCKET™ TT-S OPTIMIZES KNEE FLEXION WITHOUT SACRIFICING PERFORMANCE OR AESTHETICS, PROVIDING NEW LEVELS OF FORCE TRANSITIONS THAT DYNAMICALLY CONFORM WITH THE PATIENT'S EVERY STRIDE. THE SOCKET OFFERS SELECT AREAS OF FLEX AND RIGIDITY FOR UNPARALLELED ALL-DAY COMFORT.

1.

Designed for your life

The sleek Infinite Socket™ TT-S features a natural profile while wearing pants, yet delivers universal adjustability for both distal/ proximal volume management while providing posterior relief.

2.

Conforming to your every step

The human inspired design comprises a mixture of rigid and flexible areas that relieve bony prominences and reduces the risk of skin ulcerations. The interface delivers product durability, dynamic flex, and vibration dampening technology to compliment your every step.



3.

Unparalleled knee flexion

The adjustable posterior firmwear supports the calf, yet encourages superior knee flexion. A single release of the BOA® dial provides sitting relief for maximum comfort. Bend, kneel and squat without irritation or pinching.

4.

Adjustable without the bulk and weight

Experience exceptional control of your socket fit in a slim, lightweight design. Now you can have all the benefits of a daily adjustable socket without the associated bulk.



5.

Custom to you

Our unique fabrication method differs based on the patient's unique profile to ensure optimal socket strength and fit that compliments your lifestyle.

The socket is accompanied with a flush textile cosmetic cover so you can define your own look.

6.

Rapid donning and doffing

The adjustable BOA® dial, and posterior firmwear allow rapid donning and doffing to efficiently get the fit you want, when you want it.

7.

Built to last

Tested to ISO 10328 with an increased weight limit of 320 lbs.

"I FEEL LIKE I HAVE MORE CONTROL AND SUPPORT THAN EVER BEFORE. FOR THE FIRST TIME EVER, I GET THE SENSE THAT MY ORGANIC LEG IS CONTINUING MORE SEAMLESSLY DOWN INTO MY INORGANIC FOOT." - JOHN, F.



WE MAKE ABILITY



lim_innovations



liminnovations



liminnovations

424 Ninth Street
San Francisco, CA 94103

Phone 1-844-888-8546
info@liminnovations.com

www.liminnovations.com

The Infinite Socket™ TT-S is the prosthetic industry's thinnest and lightest dual-wall adjustable socket.

The Infinite Socket™ TT-S is ideal for patients with daily volume changes, discomfort, as the design mimics human anatomy to provide unparalleled function, adjustability and comfort.



INFINITE
SOCKET™ TT-S

THE MOST
COMFORTABLE,
ADJUSTABLE, CUSTOM
SOCKET FOR
BELOW-KNEE AMPUTEES.